

**Bhakta Kavi Narsinh Mehta
University**

**BSc Semester V
Chemistry Practicals
With effect from June-2018
[CBCS]**

Instructions to Examiners

CHEMISTRY PRACTICALS: SEMESTER V

[Total: 150 Marks]

- **Internal Evaluation [CCA]: 45 Marks**
- **External Evaluation [SEE]: 105 Marks**

CCA (Internal Evaluation): Total 45 Marks

Equal weightage to all Sections i.e.

- Organic Chemistry: 15 Marks
- Inorganic Chemistry : 15 Marks
- Physical Chemistry : 15 Marks

SEE (External Evaluation): Total 105 Marks

- Three Practical Exercises are to be given one from each section-
Organic, Inorganic and Physicochemical Exercises
 - Marks Distribution for Practicals Examination :
 - Organic Separation : 40 Marks
 - Inorganic Volumetric Estimation : 30 Marks
 - Physicochemical Exercise : 35 Marks
- Total : 105 Marks**

Note:

- **Internal Evaluation: [45 Marks] by the respective College**
- **External Evaluation: [105 Marks] by the Examiners appointed by BKNMU University for the Practical Examination**

TYBSc Semester V Chemistry Practical Examination Important Guidelines to the Examiners

Dear Colleagues

Welcome to the panel of Examiners!

Most of you are experienced and knowledgeable, hence aware and familiar with the examination related responsibilities. However having been assigned the responsibility of preparing the guidelines for the forthcoming BSc Semester V Chemistry Practical Examination 2018-19, this a humble effort to formulate certain guidelines to ensure uniform evaluation at all centres and a smooth conduct of the examination.

The pages following the instructions contain:

- Specific Requirements for each Exercise
- Marking Scheme for different Exercises

The Question Papers (Slip) for the Practical Examination and formats A (Organic), B (Inorganic), C (Physico-chemical), D (Viva) and E (Consolidated Marksheet) are enclosed in a separate booklet. The booklet also contains the format for detailed batch wise and group wise schedule of the examination to be displayed on the notice board at the examination centre.

In case the copies are not sufficient, the Examiners shall request the Lab Supervisor / Principal to make arrangement for photocopying the same.

General Instructions:

- The Examiners appointed for the BSc Semester V Chemistry Practicals must report to the Principal of the concerned College at the Examination centre at least one hour before the commencement of the Examination.
- The assessment of the practical answer sheets is to be kept confidential. The examination work shall not be carried out in presence of any other member.
- The Examiner named first is the Senior Examiner
- It shall be the responsibility of the Senior Examiner to ensure smooth conduct of the examination. The other Examiners are equally responsible and shall extend wholehearted cooperation to the Senior Examiner in the smooth conduct of the examination
- The detailed batch wise and group wise schedule of the examination at the allotted centre shall be sent to the centre by Senior Examiner at least one day prior to the

commencement of the examination. For the benefit of the candidates, the schedule shall be displayed on the notice board at the centre a day before the commencement of the examination,

- For the examination purpose the candidates shall be divided into two groups 'A Group' and 'B Group'.
- The Senior Examiner shall distribute the work to his/her colleague and guide him/her where necessary.
- The Senior Examiner shall inform the Lab Supervisor to keep the answer sheets ready for the examination. The Seat Numbers shall be written in blue ballpoint pen and the Table Number in red ballpoint pen.
- The Senior Examiner shall collect the certificate from the Lab Supervisor, indicating standardization of the laboratory equipments, before the commencement of the examination.
- On the day of the completion of the examination the senior examiner shall enter the mark in the mark sheet in Blue/ Black, seal the mark sheet and forward the same to the University immediately after the examination.
- The details in the mark sheet should be filled in English only.
- The attendance report should be filled in duplicate. One copy of the attendance report should be put in the cover along with the mark sheet.
- Another copy of the attendance report and the key forms (A-E) should be packed with the answer sheets and sent, so as to reach the University before the meeting for settling the marks.
- For authenticity, the examiners shall verify the fee receipt and ID cards of the candidates.

Some Important Instructions:

- The work distribution between the two Examiners is as follows:

Examiner 1: Physico-chemical + Inorganic Volumetric

Examiner 2: Organic Separation + Viva

- Marks Weightage

Organic Separation : 35 Marks

Inorganic Volumetric	: 30 Marks
Physico-chemical	: 30 Marks
Viva (Organic +Physico-chemical)	: 10 Marks
Total Marks	: 105 Marks

▪ Practical Schedule

DAY	Session	'A' GROUP	'B' GROUP
FIRST DAY	10.00 AM, Session-1	Organic Separation	Physico Chemical
	Session-2	Physico Chemical	Organic Separation
SECOND DAY	Session-3	Volumetric Analysis & Viva	

- The number of exercises to be performed by each candidate is 3. The Examination shall be conducted in 3 Session each of 3 hours.
- The Organic Mixture given for Organic Separation should not be repeated in a batch i.e. no two students of a batch should be given the same organic mixture.
- For Inorganic Volumetric Estimation, three expert readings are to be considered. The volume given to the expert should differ from that given to the candidates.
- In case of Physico-Chemical Exercise no two students of a batch should be given the same Exercise
- The allotment of the exercise shall be as per the lucky draw system wherein the candidate shall be asked to pick up a chit indicating the exercise number.
- The Examiners are requested to go through the Question Paper Slip before assigning the Exercise. In case of any typing mistake the Examiner shall make the necessary correction before assigning the Exercise.

BSc Semester V
[CBCS]
Chemistry Practicals

Marking Scheme

Bhakta Kavi Narsinh Mehta University
BSc Semester V Chemistry Practical Examination

Organic Separation – Exercise No. 1

Note:

The organic mixture should not be repeated in a batch i.e. no two students of a batch should be given the same organic mixture.

General Instructions:

- First instruct the students to determine the type of given mixture
- After the determination of the type of mixture, instruct the students to write the method to be used for separation
- The Examiner should check and sign the separation method written by the students
- The student may then be instructed to proceed for separation of the organic mixture

In case the candidate has determine a totally wrong type for the mixture or inspite of the efforts the candidate is not in a position to determine the type, the Examiner may guide the candidate accordingly. [In such case the marks for the determination of the type of mixture may be reduced depending upon the guidance provided by the Examiner.]

Bhakta Kavi Narsinh Mehta University
BSc Semester V Chemistry Practical Examination

Organic Separation – Exercise No. 1

Marks Distribution	Total Marks: 35
1. Type of Mixture and Separation (Note: For Solid + Liquid or Liquid + Liquid mixture, marks should be assigned if the type is determined by physical methods. However the candidate should be guided to determine the nature of each component after the separation.)	7 Marks (4+3)
2. Preliminary Tests (3 marks for each component)	6 Marks
3. Detection of Elements (3 marks for each component)	6 Marks
4. Functional Group Test (3 marks for each component)	6 Marks
5. MP/BP (1 marks for each component)	2 Marks
6. Identification of the compound (1 marks for each component)	2 Marks
7. Confirmatory Tests (2 marks for each component)	4 Marks
8. Systematic Working	2 Marks
Total Marks	35 Marks

Note:

- **For Organic Mixture refer to Appendix - I (List of Organic Mixtures)**
- **For entering the marks use Form A (Organic Separation)**

APPENDIX – I
LIST OF ORGANIC MIXTURES

Type A Acid + Phenol (A+P)

- | | | |
|---------------------|---|---------------------|
| 1. Salicylic acid | + | β – Naphthol |
| 2. Benzoic acid | + | β – Naphthol |
| 3. Cinnamic acid | + | α - Naphthol |
| 4. Phthalic acid | + | α - Naphthol |
| 5. Anthranilic acid | + | β – Naphthol |
| 6. Sulphanilic acid | + | α - Naphthol |

Type B Acid + Base (A+B)

- | | | |
|---------------------|---|--------------------------|
| 7. Anthranilic acid | + | m-nitroaniline |
| 8. Sulphanilic acid | + | p-toludine |
| 9. Salicylic acid | + | diphenylamine |
| 10. Cinnamic acid | + | p-nitroaniline |
| 11. Phthalic acid | + | α - Naphthylamine |
| 12. Benzoic acid | + | p-toludine |

Type C Acid + Neutral (A+N)

- | | | |
|----------------------|---|------------------|
| 13. Benzoic acid | + | acetanilide |
| 14. Salicylic acid | + | Benzamide |
| 15. Cinnamic acid | + | Anthracene |
| 16. Phthalic acid | + | Naphthalene |
| 17. Sulphanilic acid | + | m-dinitrobenzene |
| 18. Anthranilic acid | + | Acetanilide |

Type D Phenol + Base (P+B)

- | | | |
|-------------------------|---|--------------------------|
| 19. β – Naphthol | + | m-nitroaniline |
| 20. α - Naphthol | + | p-nitroaniline |
| 21. α - Naphthol | + | p-toludine |
| 22. β – Naphthol | + | Diphenylamine |
| 23. α - Naphthol | + | α - Naphthylamine |

Type E Phenol + Neutral (P+N)

24. β - Naphthol	+	naphthalene
25. β - Naphthol	+	acetanilide
26. α - Naphthol	+	Benzamide
27. α - Naphthol	+	Anthracene
28. β - Naphthol	+	m-dinitrobenzene

Type F Base + Neutral (B+N)

29. p-toluidine	+	Anthracene
30. m-nitroaniline	+	Acetanilide
31. p-nitroaniline	+	naphthalene
32. m-nitroaniline	+	m-dinitrobenzene
33. p-toluidine	+	Benzamide

Type G Neutral + Neutral (N+N)

34. Urea	+	Benzamide
35. Urea	+	Acetanilide

Type H. Liquid + Liquid

36. Chloroform	+	Chlorobenzene (N+N)
37. Carbon tetrachloride	+	Bromobenzene (N+N)
38. Acetone	+	Bromobenzene (N+N)
39. Ethyl acetate	+	Nitrobenzene (N+N)
40. Carbon tetrachloride	+	Aniline (N+B)
41. o- Cresol	+	Toluene (P+N)
42. Ethylacetate	+	Aniline (B+N)
43. Acetone	+	Aniline (B+N)
44. Ethanol	+	Nitrobenzene (N+N)

Type I Solid + liquid

45. Salicylic acid	+	Ethylacetate (A+N)
46. Salicylic acid	+	Acetone (A+N)
47. Benzoic acid	+	Ethanol (A+N)
48. Naphthalene	+	Acetone (N+N)

- | | | |
|----------------------|---|----------------------------|
| 49. m-dinitrobenzene | + | carbon tetrachloride (N+N) |
| 50. Acetanilide | + | Acetone (N+N) |
| 51. m-dinitrobenzene | + | Acetone (N+N) |

Note

- Everyday one [liquid + solid] mixture and two [liquid + liquid] mixture or two [liquid + solid] mixture and one [liquid + liquid] mixture should be given to the candidate
- Do not repeat any mixture on the same day

Bhakta Kavi Narsinh Mehta University
BSc Semester V Chemistry Practical Examination

Inorganic Volumetric Analysis

Exercise No 2 -10 [30 Marks]

For volumetric exercise all the standard solutions are to be prepared by the students

Ex. No. 2 Estimation of Cu^{+2} and $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ in the given $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ using 0.05M $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ solution

Preparation of the Solution: Dissolve 125 gms pure crystalline $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ in 1 litre solution

Chemical Requirement:

0.05 N $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ solution to be prepared by the students

10% KI solution

2N Na_2CO_3

2N acetic acid

Fresh starch solution

Ex. No. 3. Estimation of As^{+3} and As_2O_3 in the given As_2O_3 using 0.05M $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ solution.

Preparation of the Solution: Dissolve 32 gms sodium arsenite Na_3AsO_3 in 1 litre solution

Chemical Requirement:

0.05 N $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ solution to be prepared by the students

0.05N I_2 solution

solid NaHCO_3 , fresh starch solution.

Ex. No 4(a): Estimation of the amount of Ni^{+2} in the given $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ solution using 0.01M EDTA solution. [Direct Titration]

Preparation of the Solution: Dissolve 28gms pure crystalline $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ in 1 litre solution

Chemical Requirement:

0.01M EDTA disodium salt solution to be prepared by the students

Buffer solution: Mixture of equal volume of 1M NH_4Cl solution and 1M NH_3 solution.

Murexide indicator (solid mixture from or freshly decanted solution)

Ex. No 4(b): Estimation of the amount of Ni^{+2} in the given $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ solution using 0.01M EDTA solution. [Direct Titration]

Preparation of the Solution: Dissolve 28gms pure crystalline $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ in 1 litre solution

Chemical Requirement:

0.01M EDTA disodium salt solution to be prepared by the students

Bromopyrogallol Red indicator

Buffer solution: Mixture of equal volume of 1M NH_4Cl solution and 1M NH_3 solution.

Ex. No 4(c): Estimation of the amount of Ni^{+2} in the given $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ solution using 0.01M EDTA solution. . [Back Titration]

Preparation of the Solution: Dissolve 28gms pure crystalline $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ in 1 litre solution

Chemical Requirement:

0.01M EDTA disodium salt solution to be prepared by the students

Eriochrome Black T indicator

0.01 M MgSO_4 solution

Buffer solution (pH 10): 17.5 gm NH_4Cl + 142 ml NH_3 solution make up the volume to 250 ml with distilled water

Ex. No. 5 Estimation of the amount of Mg^{+2} and Pb^{+2} in the given solution containing a mixture of Mg^{+2} and Pb^{+2} using 0.01 M EDTA solution

Preparation of the Solution: Dissolve 16 gms Lead nitrate $\text{Pb}(\text{NO}_3)_2$ and 80 gms Magnesium nitrate $\text{Mg}(\text{NO}_3)_2$ in 1 litre solution

Chemical Requirement:

0.01M EDTA disodium salt solution to be prepared by the students

Buffer solution (pH 10): 17.5 gm NH_4Cl + 142 ml NH_3 solution make up the volume to 250 ml with distilled water

10% hexamine solution or solid powder, buffer solution (pH = 10).

Indicators : Eriochrome Black-T and Xylenol Orange

Ex. No. 6 Estimation of the amount of Ca^{+2} and Zn^{+2} in the given solution containing a mixture of Ca^{+2} and Zn^{+2} using 0.01 M EDTA solution

Preparation of the Solution: Dissolve 5.5 gms CaCl_2 and 7.0 gms ZnCl_2 in 1 litre solution

Chemical Requirement:

0.01M EDTA disodium salt solution to be prepared by the students

Buffer solution (pH 10): 17.5 gm NH_4Cl + 142 ml NH_3 solution make up the volume to 250 ml with distilled water

10% hexamine solution or solid powder, buffer solution (pH = 10).

Indicators: Eriochrome Black-T and Xylenol Orange

Ex. No. 7 Determination of the amount of Fe^{+3} & Cr^{+3} in the given solution containing a mixture of $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ and $\text{CrCl}_3 \cdot 6\text{H}_2\text{O}$ using 0.01 M $\text{Pb}(\text{NO}_3)_2$ solution

Preparation of the solution: Dissolve 5.5 gms $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ and 5.5 gms $\text{CrCl}_3 \cdot 6\text{H}_2\text{O}$ in 1 litre solution

Chemical Requirements

0.01M EDTA disodium salt solution to be prepared by the students

Indicator Xylenol Orange

10% hexamine solution or solid powder,

0.01 M $\text{Pb}(\text{NO}_3)_2$ solution to be prepared by the students

Ex. No. 8 Determination of the amount of NO_2^{-1} in the given NaNO_2 or KNO_2 solution by reduction method using 0.1 N KMnO_4 solution.

Preparation of the Solution: Dissolve 35 gms NaNO_2 or 42 gms KNO_2 in 1 litre solution

Chemical Requirement:

0.1N KMnO_4 solution to be prepared by the students

2N H_2SO_4

Ex. No. 9 To determine the amount of chloride in the given sample of water using 0.02N AgNO_3 solution

Preparation of the Solution: Sample of water

Chemical Requirement:

(Approx) 0.02N AgNO_3 solution

0.02 N NaCl

K_2CrO_4 indicator solution

Ex. No. 10 To determine the purity of NaHCO_3 in the given sample of antacid

Chemical Requirements:

Sample of antacid

0.05 N HCl

methyl orange indicator

Note: For Inorganic Volumetric Analysis the Volume of the original solution given to the candidates should be 20, 22, 24, 26, 28 and 30 ml

Inorganic Volumetric Analysis

Use Form B [Inorganic Volumetric Analysis] to enter the Marks

Marks Distribution Scheme for Exercise Nos. 2-10

Total Marks: 30

Error in the reading	One Part	Two Parts (for each part)
Difference of 0.1 ml	20 Marks	10 + 10 Marks
Difference of 0.2 ml	18 Marks	9 + 9 Marks
Difference of 0.3 ml	16 Marks	8 + 8 Marks
Difference of 0.4 ml	14 Marks	7 + 7 Marks
Difference of 0.5 ml	12 Marks	6 + 6 Marks
Difference of 0.6 ml	10 Marks	5 + 5 Marks
Difference of 0.7 ml	8 Marks	4 + 4 Marks
Difference of 0.8 ml	6 Marks	3 + 3 Marks
Difference of 0.9 ml	4 Marks	2 + 2 Marks
Difference of 1.0 ml	2 Marks	1 + 1 Marks
Correct Calculations	6 Marks	6 Marks
Systematic Working	4 Marks	4 Marks
Total Marks	30 Marks	30 Marks

Bhakta Kavi Narsinh Mehta University
BSc Semester V Chemistry Practical Examination

Physico-Chemical Exercises
Exercise No 11 -24 [30 Marks]

Exercise No. 11 Conductometry

Aim: Determine the cell constant of the given conductivity cell by using 0.1N KCl solution and also determine normality of the given HCl solution by conductometric titration

Requirements:

Solid KCl

0.5N NaOH solution

0.1 N HCl solution (as unknown concentration)

Marks Distribution

Arrangement + Systematic Working	6 Marks
Correct Cell Constant	6 Marks
Neatness –Observation table (at least 8 readings)	6 Marks
Correct Graph	6 Marks
Correct Normality	6 Marks
Total Marks	30 Marks

Exercise No 12 Conductometry

Aim: Determine the cell constant of the conductivity cell by using 0.1N KCl solution and determine the concentration of each component in the given mixture of HCl + CH₃COOH in terms of normality by conductometric titration.

Requirements:

0.1 N KCl

0.5 N NaOH

0.05 N HCl (as unknown concentration)

0.05 N CH₃COOH (as unknown concentration)

Marks Distribution

Arrangement + Systematic Working	6 Marks
Correct Cell Constant	6 Marks

Neatness –Observation table (at least 8 readings)	6 Marks
Correct Graph	6 Marks
Correct Normality	6 Marks
Total Marks	30 Marks

Exercise No. 13 Conductometry

Aim: Determine the Cell constant of the given cell and determine the normality of the given CH_3COOH by conductometric titration

Requirements:

Solid KCl solution

0.5 N NaOH

0.1 N CH_3COOH (as unknown concentration)

Marks Distribution

Arrangement + Systematic Working	6 Marks
Correct Cell Constant	6 Marks
Neatness –Observation table (at least 8 readings)	6 Marks
Correct Graph	6 Marks
Correct Normality	6 Marks
Total Marks	30 Marks

Exercise No. 14 Conductometry

Aim: Determine the Cell constant of the given conductivity cell and also determine normality of the given Ni^{+2} solution by conductometric titration.

Requirements:

0.1 M EDTA

0.1 M NiSO_4 (as unknown concentration)

Marks Distribution

Arrangement + Systematic Working	6 Marks
Correct Cell Constant	6 Marks
Neatness –Observation table (at least 8 readings)	6 Marks
Correct Graph	6 Marks
Correct Normality	6 Marks
Total Marks	30 Marks

Exercise No. 15 Conductometry

Aim: To determine the cell constant of the conductivity cell and the normality of 'X' N AgNO₃ using 0.5 N NaCl by conductometric titration

Requirements:

Solid KCl

0.5 N NaCl

0.5 N AgNO₃ (as unknown concentration)

Marks Distribution

Arrangement + Systematic Working	6 Marks
Correct Cell Constant	6 Marks
Neatness –Observation table (at least 8 readings)	6 Marks
Correct Graph	6 Marks
Correct Normality	6 Marks
Total Marks	30 Marks

Exercise No. 16 Thermodynamics

Aim: To calculate entropy of vapourisation (ΔS_v) of a given liquid by Kinetic approach i.e. from the graph of log (1/time) against (1/temperature).

Requirements

Benzene / CHCl₃ / CCl₄ / n-hexane

Marks Distribution

Systematic Working	4 Marks
Presentation & Observation	8 Marks
Graph	6 Marks
Calculations	8 Marks
Unit	2 Marks
Correct Result	2 Marks
Total Marks	30 Marks

Exercise No. 17 Refractometry

Aim: Determine the refractive index of given liquid A, B, C, D. and thus calculate the specific refractive index and molecular refractive index

Requirements:

Refractometer

Any four of the following liquids: Benzene, Toluene, Xylene, n-propyl alcohol, n-

butyl alcohol, chloroform, carbon tetrachloride, chlorobenzene

Marks Distribution

Systematic handling of Refractometer	4 Marks
Correct specific density	8 Marks
Correct refractive index	10 marks
Correct specific refractive index	4 marks
Correct specific molecular refractive index	4 Marks
Total Marks	30 Marks

Exercise No. 18 Refractometry

Aim: To prepare glycerine solution with a concentration of 10%, 5% & 2.5 % from the given glycerine solution and Determine the refractive index of given liquid of water and glycerol and thus calculate the specific refractive index and molecular refractive index.

Requirements:

Refractometer

Glycerine

Marks Distribution

Systematic handling of Refractometer	2 Marks
Correct specific density of all liquids	6 Marks
Correct refractive index of all liquids	7 Marks
Correct specific refractive index of all liquids	6 Marks
Correct specific molecular refractive index	3 Marks
Graph	3 Marks
% composition of the unknown	3 Marks
Total Marks	30 Marks

Exercise No. 19 Viscosity

Aim: Find relative and absolute viscosity of given pure liquids A, B, C & D by Ostwald's viscometer

Requirement:

Ostwald's Viscometer

Specific gravity bottle

Benzene

Toluene

Carbon tetrachloride

Chloroform

Stopwatch

Marks Distribution

Systematic handling of Viscometer & Specific gravity bottle	3 Marks
Correct density of all the four liquids	6 Marks
Readings for water & all the liquids	8 Marks
Correct viscosity of all the four liquids	3 Marks
Correct relative & absolute viscosity	6 Marks
Systematic Tabulation of Readings	4 Marks
Total Marks	30 Marks

Exercise No. 20 Viscosity

Aim: To prepare glycerine solution with a concentration of 10%, 5% & 2.5 % from the given glycerine solution and determine relative and absolute viscosity of given glycerine solution and unknown concentration by Ostwald's viscometer

Requirement:

Viscometer

Specific gravity bottle

Glycerine

Stopwatch

Marks Distribution

Systematic handling of Viscometer & Specific gravity bottle	5 Marks
Correct relative & absolute viscosity	7 Marks
Correct viscosity of all the concentrations	7 Marks
Graph	6 Marks
Concentration of the unknown solution	5 Marks
Total Marks	30 Marks

Exercise No 21. Colorimetry

Aim: Determine the amount of Ni^{+2} in the given solution by colorimetric method.

Requirements:

Ni^{+2} solution: 0.500 gm Nickel sulphate or 0.673 gms nickel ammonium sulphate in 1 litre solution

Saturated Br₂ solution,
1 % alcoholic solution of DMG (Dimethylglyoxine)
1: 1 NH₃ solution
Green Filter

Note:

In a 100 ml measuring flask give the candidate 5 to 10 ml Ni⁺² solution. As an unknown give 3, 5 or 7 ml.

Marks Distribution

Arrangement + Systematic Working	4 Marks
Preparation of dilute solutions	4 Marks
Correct OD and T %	10 Marks
Correct Graph (Two graphs)	6 Marks
Correct Concentration	6 Marks
Total Marks	30 Marks

Exercise No. 22. Colorimetry

Aim: To determine the amount of Fe⁺³ in the given solution by colorimetric method.

Requirement:

Ferric Ammonium Sulphate solution (0.3mg Fe /ml) in a 100 ml flask,
20% KCNS solution or 50 % NH₄CNS
2N HCl

Preparation of Ferric Ammonium Sulphate solution: Dissolve 2.6 gms of Ferric Ammonium Sulphate in 10ml Fe free dil HCl add distilled water and make up the volume to 1 litre (i.e. 0.3mg Fe /ml)

Preparation of 20% KCNS solution: 20 gms KCNS in 100 ml solution

Preparation of 50 % NH₄CNS solution: 50 gms NH₄CNS in 100 ml solution

Note:

In a 100 ml measuring flask give the candidate 5 to 10 ml Fe⁺³ solution. As an unknown give 3, 5 or 7 ml.

Marks Distribution

Arrangement + Systematic Working	4 Marks
Preparation of dilute solutions	4 Marks
Correct OD and T %	10 Marks
Correct Graph (Two graphs)	6 Marks

Correct Concentration	6 Marks
Total Marks	30 Marks

Exercise No. 23 Polarimeter

Aim : To determine specific rotation of the given dextrose solution of three different concentration (10%, 5%, 2.5%). Plot the graph of specific rotation against concentration and hence determine the unknown concentration.

Requirements:

Glucose, Tartaric acid or Sucrose

Marks Distribution

Arrangement + Systematic Working	4 Marks
Observations for angle of rotation & Neatness	8 Marks
Calculation of Specific Rotation	6 Marks
Graphs	6 Marks
Correct Readings	4 Marks
Correct Unknown concentration	2 Marks
Total Marks	30 Marks

Exercise No. 24 Polarimeter

Aim: To study the inversion rate of sugar in presence of 1N HCl and hence determine the rate constant for the inversion of cane sugar.

Requirements

1. Cane sugar
2. 1N HCl

Marks Distribution

Arrangement + Systematic Working	6 Marks
Observations for angle of rotation & neatness	8 Marks
Calculation of Rate constant	4 Marks
Graphs	6 Marks
Determination of order of reaction	4 Marks
Result	2 Marks
Total Marks	30 Marks

**Bhakta Kavi Narsinh Mehta
University**

**BSc Semester V
Chemistry Practicals
With effect from June-2018
[CBCS]**

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination

ORGANIC SEPARATION

Time: 3.00 Hrs

Marks: [35+5]

Exercise No. 1
Organic Separation

In the container bearing the **Capsule No** you are given a mixture of two organic compounds with distinct chemical and physical properties. The given mixture belongs to one of the types given below.

By preliminary observations/ tests decide the type of the mixture.

Get the signature of the Examiner before proceeding to separate the mixture.

Separate the two components of the mixture and identify each components of the mixture.

Perform all the tests and identify each component of the mixture

Keep both the purified components for inspection by the Examiners

Type of Mixture:.....

Final result:

No.	Nature	Elements	Functional Group	MP/BP	Name	Structural Formula
1						
2						

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
INORGANIC VOLUMETRIC ANALYSIS

Time: 3.00 Hrs

Marks: 30

Iodometry / Iodimetry

Exercise No. 2

Estimation of Cu

Aim: To determine by iodometric method the amount of Cu^{+2} and $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ in the given $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ using 0.05N $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ solution.

Requirements: Solid $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$, 10% KI, 2N Na_2CO_3 solution, 2N acetic acid, fresh starch solution.

Procedure:

Step –I: Preparation of 0.05 M $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ solution


Weigh 3.1025 gm $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ in a previously weighed watch glass and dissolve in distilled water making up the volume to 250 ml in a measuring flask. Fill this 0.05 M $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ solution in a burette.

Step –II: Estimation of Cu^{+2} and $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$

Dilute, with distilled water, the given acidic copper sulphate solution to 250 ml in a measuring flask. Take 25 ml of the diluted solution in a conical flask and add dilute Na_2CO_3 or NH_4OH drop wise with constant stirring. (To remove mineral acid) Continue adding 2 N Na_2CO_3 till complete precipitation (i.e. till the ppts no longer dissolve on stirring)

Now add 2N CH_3COOH to this solution till ppts completely dissolve. The solution turns clear blue. Add 20 ml 10% KI solution and immediately titrate the liberated iodine against $\text{Na}_2\text{S}_2\text{O}_3$ solution filled in the burette; when the solution in the flask turns light yellow add approximately 1 ml fresh starch solution (the solution turns blue). Continue adding $\text{Na}_2\text{S}_2\text{O}_3$ solution from the burette till the solution turns white (stable). The white ppt obtained is Cu_2I_2 . Note the burette reading.

Calculations:



Result

1. 25ml dil. $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ solution requires ml. 0.05 M $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ solution

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

2. The amount of Cu^{+2} in the given solution = g
3. Amount of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ in the solution = g

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
INORGANIC VOLUMETRIC ANALYSIS

Time: 3.00 Hrs

Marks: 30

Iodometry / Iodimetry

Exercise No. 3

Estimation of As

Aim: To determine by iodimetric method the amount of As^{+3} and As_2O_3 in the given As_2O_3 using 0.05M $Na_2S_2O_3 \cdot 5H_2O$ solution.

Requirements: Solid $Na_2S_2O_3 \cdot 5H_2O$, 0.05N I_2 solution, solid $NaHCO_3$, fresh starch solution.

Procedure:

Preparation of 0.05 M $Na_2S_2O_3 \cdot 5H_2O$ solution

Weigh 3.1025 gm $Na_2S_2O_3 \cdot 5H_2O$ in a previously weighed watch glass and dissolve in distilled water making up the volume to 250 ml. Fill this 0.05 M $Na_2S_2O_3 \cdot 5H_2O$ solution in a burette.

Part-I: Standardization of Iodine solution

Take 25 ml of the given iodine solution and titrate it against $Na_2S_2O_3 \cdot 5H_2O$ solution filled in the burette. When the solution turns yellow add approximately 1 ml starch solution. The solution turns blue. Add $Na_2S_2O_3 \cdot 5H_2O$ solution from the burette with constant stirring till the solution turns colourless. Note the burette reading.

Part –II: Estimation of As_2O_3

Note: Arsenic salts are highly poisonous; hence instead of using a pipette, use a burette for taking Arsenic. Avoid contact of arsenic salts if skin is cut or bruised.

Dilute with distilled water the given As_2O_3 solution to 250 ml in a measuring flask. Fill this diluted As_2O_3 solution in the burette.

Take 25 ml I_2 solution in a conical flask; add approximately 50ml (2 test tubes) distilled water and approximately 3 gms pure $NaHCO_3$ (solid). Shake the solution to dissolve the solid $NaHCO_3$. Titrate this solution against As_2O_3 solution filled in the burette. When the solution in the flask turns light yellow add approximately 1 ml fresh starch solution (the solution turns blue). Continue adding As_2O_3 solution from the burette till the solution turns colourless. Note the burette reading.

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

Calculation: Estimation of Arsenic

$$1000 \text{ ml } 1\text{N } \text{I}_2 \text{ solution} \equiv 49.45 \text{ g } \text{As}_2\text{O}_3 \equiv 37.45 \text{ g } \text{As}^{+3}$$

Result

1. Normality of I_2 solution = N
2. Volume of 'X' N I_2 solution required for 25 ml diluted As_2O_3 solution is ml
3. Amount of As_2O_3 in the given solution is g
4. Amount of As^{+3} in the given solution is g

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
INORGANIC VOLUMETRIC ANALYSIS

Time: 3.00 Hrs

Marks: 30

Complexometric

Exercise No. 4(a)

Estimation of Ni

[Direct Titration Method using Murexide as an indicator]

Aim: To determine the amount of Ni^{2+} in the given $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ solution using 0.01 M EDTA solution.

Requirements: EDTA disodium salt (solid), 1M NH_4Cl solution, liquor ammonia, distilled water, Murexide indicator (solid mixture from or freshly decanted solution).

Procedure:

Step –I: Preparation of 0.01 M EDTA disodium salt solution

Dissolve 0.9306 g EDTA salt in distilled water and make up the volume to 250 ml in a measuring flask. Fill this solution in the burette.

Step –II: Estimation of Ni^{2+} [Using Murexide as an indicator]

Dilute with distilled water the given Ni^{2+} solution to 250ml in a measuring flask. Take 25 ml of the diluted Ni^{2+} solution, add 150 ml distilled water and add about 40-50 mg (approximately 1 pinch) murexide indicator (indicator + solid KNO_3 mixture) and shake vigorously. Now add 10 ml Buffer solution (Mixture of 1 M NH_4Cl solution and 1 M Ammonia) till distinct yellow colour is obtained. Titrate the solution against 0.01M EDTA disodium salt solution with constant stirring. At the end point the yellow colour of the solution turns violet. Note the burette reading

Calculations

$$1000 \text{ ml } 1\text{M EDTA-salt} \quad \equiv \quad 58.7 \text{ g } \text{Ni}^{2+} \quad \equiv \quad 280.8 \text{ g } \text{NiSO}_4 \cdot 7\text{H}_2\text{O}$$

Result

1. Volume of 0.01M EDTA solution required for 25 ml diluted Ni^{2+} solution is ml
2. Amount of Ni^{2+} in the given solution is g
3. Amount of $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ in the given solution is g

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Inorganic Volumetric Analysis

Time: 3.00 Hrs

Marks: 30

Complexometric

Exercise No. 4(b)

Estimation of Ni

Direct Titration Method using Bromopyrogallol Red as an indicator

Aim: To determine the amount of Ni⁺² in the given NiSO₄·7H₂O solution using 0.01 M EDTA solution.

Requirements: EDTA disodium salt (solid), Buffer solution, Bromopyrogallol Red indicator, distilled water

Buffer solution: Mixture of equal volume of 1M NH₄Cl solution and 1M NH₃ solution.

Procedure:

Step –I: Preparation of 0.01 M EDTA disodium salt solution

Dissolve 0.9306 g EDTA salt in distilled water and make up the volume to 250 ml in a measuring flask. Fill the solution in the burette.

Step –II: Estimation of Ni²⁺ [Using Bromopyrogallol Red as an indicator]

Dilute with distilled water the given Ni²⁺ solution to 250 ml in a measuring flask. Take 25 ml of the diluted Ni²⁺ solution, add 150 ml distilled water and 1 ml Bromopyrogallol Red indicator. Now add 10 ml buffer solution and shake vigorously. Titrate the solution against 0.01M EDTA disodium salt solution from the burette. At the end point the blue colour of the solution turns red. Note the burette reading

Calculations

$$1000 \text{ ml } 1\text{M EDTA-salt} \quad \equiv \quad 58.7 \text{ g Ni}^{2+} \quad \equiv \quad 280.8 \text{ g NiSO}_4 \cdot 7\text{H}_2\text{O}$$

Result

1. Volume of 0.01M EDTA solution required for 25 ml diluted Ni²⁺ solution is ml
2. Amount of Ni²⁺ in the given solution is g
3. Amount of NiSO₄·7H₂O in the given solution is g

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Inorganic Volumetric Analysis

Time: 3.00 Hrs

Marks: 30

Complexometric

Exercise No. 4(c)

Estimation of Ni

[Back Titration Method using Eriochrome Black T indicator]

Aim: To determine the amount of Ni^{2+} in the given $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ solution using 0.01 M EDTA solution.

Requirements: EDTA disodium salt (solid), Buffer solution pH 10, Eriochrome Black T indicator, 0.01 M MgSO_4 solution, distilled water

Procedure

Step –I: Preparation of 0.01 M EDTA disodium salt solution

Dissolve 0.9306 g EDTA salt in distilled water and make up the volume to 250 ml

Step –II: Preparation of 0.01 M (exact) MgSO_4 solution

Dissolve 0.616 g MgSO_4 in distilled water and make up the volume to 250 ml

Step –III: Blank Reading

Fill the prepared 0.01 M MgSO_4 solution in the burette. Take 25 ml EDTA solution in a conical flask add 25 ml distilled water and 3 ml buffer solution. Add 2-3 drops Eriochrome Black T indicator and titrate against 0.01 M MgSO_4 solution from the burette. At the end point the blue coloured solution turns wine red.

Step – IV: Estimation of Ni^{2+}

Dilute with distilled water the given Ni^{2+} solution to 250ml in a measuring flask. Take 25 ml of the diluted Ni^{2+} solutions (if the solution is acidic neutralize by adding NH_4OH – then heat to remove excess NH_3 ; cool the flask under tap water). Add a definite volume (25 ml) of 0.01M EDTA solution. Now add 25 ml distilled water, 3 ml buffer solution and 2-3 drops Eriochrome Black T indicator. Shake and titrate the solution against MgSO_4 solution from the burette. At the end point the blue colour of the solution turns red. Note the burette reading

Calculation

$$1000 \text{ ml } 1\text{M EDTA-salt} \equiv 58.7 \text{ g Ni}^{2+} \equiv 280.8 \text{ g NiSO}_4 \cdot 7\text{H}_2\text{O}$$

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

Result

1. Volume of 0.01M EDTA solution required for 25 ml diluted Ni^{2+} solution is ml
2. Amount of Ni^{2+} in the given solution is g
3. Amount of $\text{NiSO}_4 \cdot 7\text{H}_2\text{O}$ in the given solution is g

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Inorganic Volumetric Analysis

Time: 3.00 Hrs

Marks: 30

Complexometric

Exercise No. 5

Estimation of Pb and Mg

Aim: To determine the amount of Mg^{+2} and Pb^{+2} in the given solution containing a mixture of Mg^{+2} and Pb^{+2} using 0.01 M EDTA solution

Requirements: EDTA– disodium salt (solid), distilled water, indicators Eriochrome Black-T and Xylenol Orange, 10% hexamine solution or solid powder, buffer solution (pH = 10).

Procedure:

Step -I: Preparation of 0.01M EDTA-disodium salt solution:

Dissolve the 0.9306 gm EDTA in distilled water & make up the total volume to 250 ml

Step- II: Dilute with distilled water the given solution containing a mixture of Mg^{+2} and Pb^{+2} to 250 ml with distilled water.

Step –III: Determination of the total amount of Pb^{+2} and Mg^{+2} in the given mixture

Take 25 ml of the diluted mixture in a conical flask. Add 50 ml distilled water and 5 ml (definite volume) of buffer solution (pH=10). Now add 3-4 drops Eriochrome Black-T as an indicator and titrate against 0.01M EDTA solution from the burette with constant stirring. At the end point the wine red coloured solution turns blue. Note the burette reading.

Step – IV: Determination of the amount of Pb^{+2} in the solution

Take 25ml diluted mixture in a conical flask and add 50ml distilled water and 3-4 drops Xylenol Orange indicator. Now add 3 ml dil. HNO_3 (till solution turns yellow). Then slowly add 5 ml 10% hexamine solution or the solid powder with constant shaking till the yellow solution turns red (pH=6). Now titrate the solution against 0.01M EDTA solution from the burette. At the end point the red solution turns distinct yellow. Note the burette reading.

Calculations

1000 ml 1M EDTA-salt \equiv 207.2 g Pb^{2+} \equiv 331.2 g $Pb(NO_3)_2$

1000 ml 1M EDTA-salt \equiv 24.3 gm Mg^{2+} \equiv 148.3 gm $Mg(NO_3)_2$

Result

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

1. Amount of Pb^{+2} in the given solution is g
2. Amount of $\text{Pb}(\text{NO}_3)_2$ in the given solution is g
3. Amount of Mg^{+2} in the given solution is g
4. Amount of $\text{Mg}(\text{NO}_3)_2$ in the given solution is g

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Inorganic Volumetric Analysis

Time: 3.00 Hrs

Marks: 30

Complexometric

Exercise No. 6

Estimation of the amount of Ca and Zn

Aim: To determine the amount of Ca^{+2} & Zn^{+2} in the given solution containing a mixture of $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ and ZnCl_2 using 0.01 M EDTA solution

Requirements: EDTA – disodium salt (solid), distilled water, indicators Eriochrome Black-T and Xylenol Orange, 10% hexamine solution or solid powder, buffer solution (pH = 10).

Procedure:

Step -I: Preparation of 0.01M EDTA-disodium salt solution:

Dissolve the 0.9306 gm EDTA in distilled water & make up the total volume to 250 ml

Step- II: Dilute with distilled water the given solution containing a mixture of Ca^{+2} and Zn^{+2} to 250 ml with distilled water.

Step –III: Determination of the total amount of Ca^{+2} and Zn^{+2} in the given mixture

Take 25 ml of the diluted mixture in a conical flask. Add 25 ml distilled water and 5 ml (definite volume) of buffer solution (pH= 10). Now add 3-4 drop Eriochrome Black-T as an indicator and titrate against 0.01M EDTA solution from the burette with constant stirring. At the end point the wine red coloured solution turns blue. Note the burette reading.

Step – IV: Determination of the amount of Zn^{+2} in the solution

Take 25ml diluted mixture in a conical flask and add 25 ml distilled water and 3-4 drops Xylenol Orange indicator (solution turns yellow). Now add 5 ml 10% hexamine solution drop wise or slowly add the solid powder with constant shaking till the yellow solution turns red (pH=6). Now titrate the solution against 0.01M EDTA solution from the burette. At the end point the red solution turns distinct yellow. Note the burette reading.

Calculations

1000 ml 1M EDTA-salt \equiv 65.4 g Zn^{+2} \equiv 136.4 g ZnCl_2

1000 ml 1M EDTA-salt \equiv 40 gm Ca^{2+} \equiv 147 gm $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$

Result

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

1. Amount of Zn^{+2} in the given solution is g
2. Amount of ZnCl_2 in the given solution is g
3. Amount of Ca^{+2} in the given solution is g
4. Amount of $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ in the given solution is g

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Inorganic Volumetric Analysis

Time: 3.00 Hrs

Marks: 30

Complexometric

Exercise No. 7

Estimation of the amount of Fe and Cr

Aim: To determine the amount of Fe^{+3} & Cr^{+3} in the given solution containing a mixture of $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ and $\text{CrCl}_3 \cdot 6\text{H}_2\text{O}$ using 0.01 M $\text{Pb}(\text{NO}_3)_2$ solution

Requirements: EDTA – disodium salt (solid), distilled water, indicators Xylenol Orange, 10% hexamine solution or solid powder, 0.01 M $\text{Pb}(\text{NO}_3)_2$.

Procedure:

Step -I: Preparation of 0.01M EDTA-disodium salt solution:

Dissolve 0.9306 gm EDTA in distilled water & make up the total volume to 250 ml.

Step – II: Preparation of 0.01M $\text{Pb}(\text{NO}_3)_2$ solution:

Dissolve the 0.828 gm $\text{Pb}(\text{NO}_3)_2$ in distilled water & make up the total volume to 250 ml.

Step- III: Dilute with distilled water the given solution containing a mixture of Fe^{+3} and Cr^{+3} to 250 ml with distilled water.

Step –III: Determination of the total amount of Fe^{+3} in the given mixture

Take 25 ml of the diluted solution of mixture containing Fe^{+3} and Cr^{+3} in a conical flask. Add 50 ml distilled water and 25 ml 0.01M EDTA solution. Now add 12-15 ml 10% hexamine solution (pH between 5-6 use pH paper) and 3-4 drops xylenol orange as an indicator (solution turns yellow) and titrate against 0.01M $\text{Pb}(\text{NO}_3)_2$ solution from the burette with constant stirring. At the end point the yellow coloured solution turns first orange then red-violet. **(Do not throw away this solution and use it for part-II)** Note the burette reading (A ml).

Step – IV: Determination of the amount of Cr^{+3} in the given mixture

To the red-violet solution obtained in part I, add 25 ml 0.01 M EDTA solution and 10 ml 1N HNO_3 solution (pH should be between 1 & 2). Boil the solution for 15 to 20 min, cool and add 25 ml distilled water. Now add 10 ml 10% hexamine solution and 3-4 drops xylenol orange as an indicator. The solution turns yellow. Now titrate the solution against 0.01M $\text{Pb}(\text{NO}_3)_2$ solution from the burette with constant stirring. At the end point the yellow solution will turn orange and then red-violet. Note the burette reading (B ml).

Calculations

1000 ml 1M EDTA-salt \equiv 55.83 g Fe^{+3} \equiv 270.3 g $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

1000 ml 1M EDTA-salt \equiv 52 g Cr⁺³ \equiv 266.5 g CrCl₃.6H₂O

Result

1. Amount 0.01M EDTA solution used for Fe⁺³ in 25 ml diluted solution of the mixture _____ (25-A) ml
2. Amount 0.01M EDTA solution used for Cr⁺³ in 25 ml diluted solution of the mixture _____(25-B) ml
3. Amount of Fe⁺³ in the given solution is g
4. Amount of Cr⁺² in the given solution is g

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Inorganic Volumetric Analysis

Time: 3.00 Hrs

Marks: 30

Redox

Exercise No. 8

Estimation of NO_2^{-1} ions

Aim: To determine, by redox titration, the amount of NO_2^{-1} ions in the given NaNO_2 solution, with the help of 0.1N KMnO_4 solution.

Requirement: Solid KMnO_4 , 2N H_2SO_4

Procedure:

Step – I : Preparation of 0.1N KMnO_4 solution.

Dissolve 0.79gm of KMnO_4 in distilled water and make the volume to 250 ml.

Step –II: Estimation of NO_2^{-1}

Dilute, with distilled water, the given NaNO_2 solution to 250ml in a measuring flask. Fill this solution in the burette.

Now take 25ml of 0.1 N KMnO_4 solution in a beaker. Add 25ml of 2N H_2SO_4 and 25ml distilled water. Heat it upto 50°C and then titrate against NaNO_2 from burette. [Precaution: During titration take care to see that the jet of the burette is dipped in the beaker used for titration] At the endpoint the light pink colour becomes colourless. Note the burette reading.

Calculation

$$1000\text{ml } 1\text{N } \text{KMnO}_4 \equiv 23 \text{ gm } \text{NO}_2^{-1} \quad \equiv \quad 34.5 \text{ gm } \text{NaNO}_2$$

Result

1.Amount of NO_2^{-1} in the given solution = gm

2.Amount of NaNO_2 in the given solution = gm

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Inorganic Volumetric Analysis

Time: 3.00 Hrs

Marks: 30

Water analysis

Exercise No. 9

Determination of Chloride

Aim: To determine the amount of chloride in the given sample of water using 'X' N AgNO₃ solution.

Requirements: 'X' N AgNO₃ solution, 0.02 N NaCl, K₂CrO₄ indicator solution

Procedure

Standardization of AgNO₃

Take 10 ml 0.02 N NaCl in a porcelain dish, add 2 drops K₂CrO₄ indicator. Titrate, with constant stirring, against 'X' N AgNO₃ solution until a slight reddish colour ppt is obtained. Note the burette reading.

Estimation of Cl¹⁻ ions in the water sample (tap water)

Take 25 ml of the water sample (tap water) in a porcelain dish, add 2 drops K₂CrO₄ indicator. Titrate, with constant stirring, against 'X' N AgNO₃ solution from the burette until a slight reddish colour ppt is obtained. Note the burette reading.

Calculations

1000 ml 1N AgNO₃ ≡ 35.5 gms Chloride

Result:

1. The Normality of the given AgNO₃ solution is N
2. Volume of AgNO₃ solution required for 25 ml sample of water is ml
3. Amount of chloride in the given sample of water is mg/litre

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Inorganic Volumetric Analysis

Time: 3.00 Hrs

Marks: 30

Exercise No. 10

Determination of Purity of NaHCO₃

Aim: To determine the purity of NaHCO₃ in the given sample of antacid by volumetric analysis.

Requirements: Sample of antacid, 0.05 N HCl, methyl orange indicator.

Procedure:

Method: Dissolve the given sample of antacid in distilled water and make up volume to 250ml. If necessary filter so as to get a clear solution.

From that pipette out 25ml solution in conical flask and add 2-3 drops methyl orange indicator and titrate against 0.05N HCl solution. At the end point color change will be from yellow to orange. (A) ml.

Calculations:

Amount of 0.05N HCl solution used for 25ml sample solution = _____(A) ml.

For 250ml diluted solution used up amount of 0.05N HCl = A × 10 = _____ (B) ml

Ask the Examiner for the weight of the sample of antacid

Practically equivalent wt. of antacid (Ep) = $\frac{\text{Wt. of antacid} \times 1000}{\text{'B' ml} \times \text{'N' of HCl}}$

Theoretically equivalent wt. of antacid (Ec) = 84 gm for 100 % pure

Practically equivalent wt. of antacid (Ep) = $\frac{Ep \times 100}{Ec}$
= 'X'
= 'X' % purity of antacid.

Result:

1. Amount of 0.05N HCl used for the neutralization of antacid= _____ (A) ml
2. Practically equivalent weight (Ep) of antacid = _____ (Ep)
3. Purity of antacid = _____ (X) %.

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 11 Conductometry

Aim: Determine the cell constant of the given conductivity cell by using 0.1N KCl solution and also determine normality of the given HCl solution by conductometric titration

Requirements: 0.5N NaOH, 'X'N HCl, solid KCl.

Procedure:

1. Prepare 100ml 0.1N KCl solution and use it to determine the cell constant of the conductivity cell.
2. Take 25ml of the given HCl solution in a beaker (or use the volume as per the capacity of the cell) and dip the conductivity cell in it. Determine the conductivity of the cell and also determine the given HCl and the conductivity after each addition of 0.5 ml NaOH, till the neutralization point is reached. Take similar readings after the neutralization point is reached.

Observation No	Volume of added 0.5N NaOH solution [V ml]	Conductance 'C' mho
1	0.0 ml	
2	0.5 ml	
3	1.0 ml	
4	1.5 ml	
5	2.0 ml	
6	2.5 ml	
7	3.0 ml	
8	3.5 ml	
9	4.0 ml	
10	4.5 ml	
11	5.0 ml	
12	5.5 ml	
13	6.0 ml	
14	6.5 ml	
15	7.0 ml	
16	7.5 ml	

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

17	8.0 ml	
18	8.5 ml	
19	9.0 ml	
20	9.5 ml	

Graph: Plot the graph of conductance against the volume of alkali added and use it to determine the point of neutralization. Calculate the normality of the given HCl solution.

- Result:**
1. Cell constant of the conductivity cell = -----
 2. Normality of HCl =

BHAKTA KAVI NARSINH MEHTA UNIVERSITY**BSc Semester V****Chemistry Practical Examination****Physico Chemical Exercise****Time: 3.00 Hrs****Marks: [30 +5]****Exercise No. 12 Conductometry**

Aim: Determine the cell constant of the conductivity cell by using 0.1N KCl solution and determine the concentration of each component in the given mixture of HCl + CH₃COOH in terms of normality by conductometric titration.

Requirements: 0.5N NaOH, XN HCl, 'Y' N CH₃COOH, 0.1 N KCl

Procedure: Take 25ml of the given 'X' N HCl + 25 ml 'Y' N CH₃COOH solution in a beaker (or use the volume as per the capacity of the cell) and dip conductivity cell in it. Determine the conductivity of the given mixture of HCl + CH₃COOH solution. Similarly take more readings by adding 0.5 ml 0.5 N NaOH solution. Take about 30 readings. Continue the titration and take the reading till it gives the conductance for the strong alkali

Observation No	Vol. of added 0.5N NaOH solution [V ml]	Conductance 'C' mho
1	0.0 ml	
2	0.5 ml	
3	1.0 ml	
4	1.5 ml	
5	2.0 ml	
6	2.5 ml	
7	3.0 ml	
8	3.5 ml	
9	4.0 ml	
10	4.5 ml	
11	5.0 ml	
12	5.5 ml	
13	6.0 ml	
14	6.5 ml	
15	7.0 ml	
16	7.5 ml	
17	8.0 ml	
18	8.5 ml	
19	9.0 ml	
20	9.5 ml	

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

Graph: Plot the graph of volume of alkali added against conductance and use it to determine the point of neutralization.

Calculate the normality of each component in the given mixture of HCl + CH₃COOH.

Result

1. Cell constant of the conductivity cell = -----
2. Normality of HCl =
3. Normality of CH₃COOH =

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 13 Conductometry

Aim: Determine the Cell constant of the given conductivity cell and determine the normality of the given CH_3COOH by conductometric titration

Requirements: 0.5N NaOH, 'X' N CH_3COOH , Prepare 0.1 N KCl solution

Procedure: Take 25ml of the given 'X' N CH_3COOH solution in a beaker (or use the volume as per the capacity of the cell) and dip conductivity cell in it. Determine the conductivity of the given X N CH_3COOH solution. Similarly take more readings by adding 0.5 ml 0.5 N NaOH solution. Take about 30 readings.

Observation No	Vol of added 0.5N NaOH solution [V ml]	Conductance 'C' mho
1	0.0 ml	
2	0.5 ml	
3	1.0 ml	
4	1.5 ml	
5	2.0 ml	
6	2.5 ml	
7	3.0 ml	
8	3.5 ml	
9	4.0 ml	
10	4.5 ml	
11	5.0 ml	
12	5.5 ml	
13	6.0 ml	
14	6.5 ml	
15	7.0 ml	
16	7.5 ml	
17	8.0 ml	
18	8.5 ml	
19	9.0 ml	
20	9.5 ml	

Graph: Plot the graph of conductance against the volume of alkali added and use it to determine the point of neutralization. Calculate the normality of the given CH_3COOH solution

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination

Result

1. Cell constant of the conductivity cell = -----
2. Normality of CH_3COOH =

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 14 Conductometry

Aim: Determine the Cell constant of the given conductivity cell and also determine normality of the given Ni^{+2} solution by conductometric titration.

Requirements: 0.1 M EDTA, 0.1 M NiSO_4

Procedure: Take 50 ml of the given 0.1 M NiSO_4 solution in a beaker (or use the volume as per the capacity of the cell) and dip conductivity cell in it. Determine the conductivity of the given NiSO_4 solution and also determine the conductivity after each addition of 0.5 ml EDTA. Take similar readings after the end point is reached.

Observation No	Vol of added 0.1 M EDTA solution [V ml]	Conductance 'C' mho
1	0.0 ml	
2	0.5 ml	
3	1.0 ml	
4	1.5 ml	
5	2.0 ml	
6	2.5 ml	
7	3.0 ml	
8	3.5 ml	
9	4.0 ml	
10	4.5 ml	
11	5.0 ml	
12	5.5 ml	
13	6.0 ml	
14	6.5 ml	
15	7.0 ml	
16	7.5 ml	
17	8.0 ml	
18	8.5 ml	
19	9.0 ml	
20	9.5 ml	

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination

Graph: Plot the graph of conductance against the volume of EDTA added and use it to determine the point of end point. Calculate the normality of the given NiSO₄ solution.

Result

1. Cell constant of the conductivity cell = -----
2. Normality of NiSO₄ = N

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 15 Conductometry

Aim: To determine the cell constant of the conductivity cell and the normality of 'X' N AgNO₃ using 0.5 N NaCl by conductometric titration

Requirements: 0.5N NaCl , 'X' N AgNO₃, Solid KCl

Method: Prepare 100 ml solution of 0.1 N KCl and use it to determine the cell constant of the conductivity cell.

Take 25 ml 'X' N AgNO₃ solution in a beaker and dip the conductivity cell in it. Titrate the solution against 0.5 N NaCl solution. Similarly take more readings by adding 0.5 ml 0.5 N NaCl solution. After the neutralization point take more reading to determine the end point accurately.

Observation Table:

Observation No.	Vol of added 0.5N NaOH solution [V ml]	Conductance 'C' mho
1	0.0 ml	
2	0.5 ml	
3	1.0 ml	
4	1.5 ml	
5	2.0 ml	
6	2.5 ml	
7	3.0 ml	
8	3.5 ml	
9	4.0 ml	
10	4.5 ml	
11	5.0 ml	
12	5.5 ml	
13	6.0 ml	
14	6.5 ml	
15	7.0 ml	
16	7.5 ml	
17	8.0 ml	

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

18	8.5 ml	
19	9.0 ml	
20	9.5 ml	

Graph: Plot the graph of conductance against the volume of NaCl added and use it to determine the point of neutralization. Calculate the normality of the given AgNO₃ solution

Result

1. Cell Constant = _____
2. Normality of AgNO₃ = _____ N

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 16 Thermodynamics

Aim: To calculate entropy of vapourisation (ΔS_v) of a given liquid by Kinetic approach i.e. from the graph of $\log(1/t)$ against $(1/\text{temperature})$.

Requirements

Benzene / CHCl_3 / CCl_4 / n-hexane

Constant temperature bath

Procedure

With the help of a micropipette, take 0.1ml of the given liquid [Benzene / CHCl_3 / CCl_4 / n-hexane] and transfer it to a clean evaporating dish kept floating in a constant temperature bath. Note the time for complete evaporation of the liquid. Repeat the process atleast thrice.

Similarly perform the experiment at four different temperatures (with a difference of 5°C)

Tabulate your readings as follows:

Temperature for Experiment = $T^\circ\text{C}$

No	Temp $^\circ\text{C}$	Temp T $^\circ\text{K}$	1/T Express in the power of 10^{-3}	Time in seconds (t)				Log 1/t	ΔH_v kcal/mol
				t_1	t_2	t_3	Average 't' sec		
1									
2									
3									
4									
5									

Note: Ask the Examiner for the boiling point of the given liquid

Graph

- ✓ Plot the graph of $\log 1/t$ against $1/T$
- ✓ From the slope determine the Heat of Vapourization and hence calculate Entropy of Vapourization

Result

1. Heat of Vapourization of the given liquid is _____ Kcal/mol deg K

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

2. Entropy of Vapourization of the given liquid is _____ cal/mol deg K

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 17 Refractometry

Aim: Determine the refractive index of given liquid A, B, C, D. and thus calculate the specific refractive index and molecular refractive index.

Procedure: Clean specific gravity bottle with distilled H₂O, alcohol, ether and dry with dryer.

Weigh given liquid A, B, C & D and water and find out specific density and absolute density. Arrange the instrument in sunlight expose the prism and clean with water, alcohol and ether in that order and then with cotton. Lastly clean the prism with the liquid to be used for the experiment. Then put 2-3 drops of that liquid on prism and close prism. Insert the liquid by glass tube or pipette through the hole on the side of prism box. Adjust the mirror in such a manner that by the reflection of light rays maximum surface of prism is brightened. Rotate the prism with the screw so half of the circle gets darken remaining becomes brighter which can be observed through telescope. If colour bands are observed than rotate the compensator on the prism to remove colour bands. Rotate the prism screw so that the limit of the bright and dark parts touches at the centre point of cross wire inside the telescope. At this time read the refractive index on the scale, note the reading upto four digits. Similarly take three readings of each liquid, determine average and calculate the specific refractive index and molecular refractive index.

Observation Table:

No.	Liquid	Specific Density	Absolute Density	Refractive Index				Specific Refractive Index	Molecular Refractive Index
				1	2	3	Average		
1	A								
2	B								
3	C								
4	D								
5	H ₂ O								

Note: Ask the Examiner for the Molecular Weight of the Liquids

Result:

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

Liquid	A	B	C	D	Water
Specific Refractive Index					
Molecular Refractive Index					

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 18 Refractometry

Aim: To prepare glycerine solution with a concentration of 10%, 5% & 2.5 % from the given glycerine solution and Determine the refractive index of given liquid of water and glycerol and thus calculate the specific refractive index and molecular refractive index.

Procedure:

Clean specific gravity bottle with distilled H₂O, alcohol, ether and dry with dryer.

Weigh 10%, 5% & 2.5 % glycerine solutions and water and determine the specific density and absolute density. Arrange the instrument in sunlight expose the prism and clean with water, alcohol and ether in that order and then with cotton. Lastly clean the prism with the liquid to be used for the experiment. Then put 2-3 drops of that liquid on prism and close prism. Insert the liquid by glass tube or pipette through the hole on the side of prism box. Adjust the mirror in such a manner that by the reflection of light rays maximum surface of prism is brightened. Rotate the prism with the screw so half of the circle gets darken remaining becomes brighter which can be observed through telescope. If colour bands are observed than rotate the compensator on the prism to remove colour bands. Rotate the prism screw so that the limit of the bright and dark parts touches the centre point of cross wire inside the telescope. At this time read the refractive index on the scale, note the reading upto four digits. Similarly take three readings of each liquid, determine average and calculate the specific refractive index and molecular refractive index.

Observation Table:

No.	Liquid	Density	Refractive Index	Specific Refractive Index	Molecular Refractive Index
1	A (Water)				
2	B (Glycerol)				
3	10 %				
4	5 %				
5	2.5 %				
6	Unknown				

Graph

Plot the graph of % composition against specific refractive index.

Determine the concentration of the unknown solution from the graph

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination

Result:

No.	Liquid	Specific Refractive Index	Molecular Refractive Index
1	A (Water)		
2	B (Glycerol)		
3	10 %		
4	5 %		
5	2.5 %		

2. Concentration of the unknown solution is

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 19 Viscosity

Aim: Find relative and absolute viscosity of given pure liquids A, B, C & D by Ostwald's viscometer

Requirement: Viscometer, Liquid A, B, C & D and distilled H₂O, Stopwatch, Specific gravity bottle

Procedure:

1. Clean specific gravity bottle respectively with chromic acid, water, alcohol and ether. Dry it with dryer. Find specific density and absolute density by weighing liquid A, B, C & D in specific gravity bottle.

2. Clean viscometer with chromic acid, water, alcohol and ether. Dry it with dryer.

Take the definite volume of liquid A [e.g. 10 ml or 15ml] in the broad part of viscometer. Suck the liquid by rubber tube in the capillary upto the upper marker. Then allow the liquid to flow up to the lower mark. Note the time in seconds for flowing the liquid down through both the marks. Use stopwatch for this. Similarly for liquid B, C, D and distilled H₂O take three readings for each. The average of all the three readings is to be used for calculation. Calculate the relative and absolute viscosity.

Observation

No	Liquid	Specific Density	Absolute Density	Time in secs			Average 't' secs	Relative Viscosity	Absolute Viscosity
				t ₁	t ₂	t ₃			
1	A								
2	B								
3	C								
4	D								
5	H ₂ O								

Result:

1) Relative Viscosity of the liquids :

A=

B=

C=.....

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination

D=

2) Absolute Viscosity of the liquid:

A=

B=.....

C=.....

D=.....

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 20 Viscosity

Aim: To prepare glycerine solution with a concentration of 10%, 5% & 2.5 % from the given glycerine solution and determine relative and absolute viscosity of given glycerine solution and unknown concentration by Ostwald's viscometer

Requirement: Viscometer, Glycerine and distilled H₂O, Stopwatch, Specific gravity bottle

Procedure:

1. Clean specific gravity bottle respectively with chromic acid, water, alcohol and ether. Dry it with dryer. Find specific density and absolute density by weighing liquid A, B, C & D in specific gravity bottle.

2. Clean viscometer with chromic acid, water, alcohol and ether. Dry it with dryer.

Take the definite volume [10 or 20 ml] of 10 % glycerine solution in the broad part of viscometer. Suck the solution a little above the upper mark by means of a rubber tube joined to the capillary. Then release the liquid to flow down to the lower mark. Use stopwatch to note the time in seconds for the liquid to flow down through both the marks. Similarly find the flowing time for 10%, 5% & 2.5 % for liquid B, C, D and distilled H₂O take three readings for each. The average of all the three readings is to be used for calculation. Calculate the relative and absolute viscosity.

Observation

No	Liquid Conc	Specific Density	Absolute Density	Time in secs			Average 't' secs	Relative Viscosity	Absolute Viscosity
				t ₁	t ₂	t ₃			
1	10 %								
2	5 %								
3	2.5 %								
4	unknown								

Graph

Plot Graph of Concentration against Absolute Viscosity and determine the concentration of the unknown from the graph

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination

Result:

1. Concentration of the unknown solution =

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 21 Colorimetry

Aim: Determine the amount of Ni^{+2} in the given solution using DMG by colorimetric method.

Requirements: Ni^{+2} solution, Saturated Br_2 solution, 1% alcoholic solution of DMG (Dimethylglyoxine), 1: 1 NH_3 solution

Procedure: Dilute with distilled water the given Ni^{+2} solution in 100ml measuring flask.

Take 2, 4, 6, 8 and 10ml of the above Ni solution in five different measuring flask each of 100 ml capacity. Add 2ml of saturated Br_2 solution, 2ml of 1: 1 NH_3 solution and 5ml DMG solution into each of the 100ml measuring flask and make up the volume to 100ml with distilled H_2O . Similarly add 2ml of saturated Br_2 solution, 2ml of 1: 1 NH_3 solution and 5ml DMG solution into the unknown solution given in a 100 ml measuring flask and make up the volume to 100ml. Shake each flask well and keep them for 10 - 15 minutes.

Standardize the calorimeter for zero and 100 positions using with distilled H_2O as standard.

[Use green filter or 480 μ wavelength] Measure the transmittance T% and optical density (O.D) of all the above prepared solution and also for the unknown solution.

Observation No.	Volume taken ml	Concentration C mg/ml	Optical Density OD	Transmittance T %
1	2.0 ml			
2	4.0 ml			
3.	6.0 ml			
4	8.0 ml			
5	10.0 ml			
6	Unknown			

Graphs

Draw the graph of

1. Concentration against T%

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

2. Concentration against OD

Find the concentration of unknown solution from the graphs.

Note: 1 ml NiSO₄.7H₂O contains 0.01mg of Ni

Result

Concentration of unknown solution:

1. From graph I concentration = _____ mg/ml

2. From graph II concentration = _____ mg/ml

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 22 Colorimetry

Aim: To determine the amount of Fe^{+3} in the given solution by colorimetric method.

Requirement: Ferric Ammonium Sulphate solution (0.3mg Fe /ml) in a 100 ml flask, 20% KCNS solution, 2N HCl, Dropping pipette/burette, 100ml measuring flask

Method: Dilute with distilled water the given Fe^{+3} solution to 100ml in measuring flask. From this diluted solution take 2, 4, 6, 8 and 10ml in five different measuring flask each of 100ml capacity. Add 1ml 2N HCl solution and 5ml 20% KCNS solution into each measuring flask and dilute it up to 100ml with distilled H_2O . Similarly add 1ml 2N HCl and 5ml 20% KCNS solution into unknown solution given in a 100 ml measuring flask and make up the volume to 100ml. Shake each flask well and keep them for 10 -15 minutes. Keep shaking each flask after few minutes.

Measure the Optical Density and Transmittance of the known and unknown solution.

[Use green filter or 480 μ wavelength]

Observation No.	Volume taken ml	Concentration C mg	Optical Density OD	Transmittance T %
1	2.0 ml			
2	4.0 ml			
3.	6.0 ml			
4	8.0 ml			
5	10.0 ml			
6	Unknown			

Graphs

Draw the following graphs

1. Graph of Transmittance against concentration
2. Graph of Optical density against concentration

Find the concentration of unknown solution from the graph and note the results

Result

Concentration of unknown solution:

1. From graph I concentration = _____ mg/ml

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

2. From graph II concentration = _____ mg/ml

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 23 Polarimeter

Aim : To determine specific rotation of the given dextrose solution of three different concentration (10%, 5%, 2.5%). Plot the graph of specific rotation against concentration and hence determine the unknown concentration.

Candidates should prepare 10%, 5% and 2.5% concentration of the solution.

No.	Concentration	Reading of Distilled water	Reading of solution	Angle of Deviation θ	Specific Rotation α
1	10 %				
2	5 %				
3	2.5 %				
4	Unknown				

Note

1. Systematic and neat entry of the results will be considered while assigning marks
2. Details of the calculations should be clearly shown

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination
Physico Chemical Exercise

Time: 3.00 Hrs

Marks: [30 +5]

Exercise No. 24 Polarimeter

Aim: To study the inversion rate of sugar in presence of 1N HCl and hence determine the rate constant for the inversion of cane sugar.

Requirements

1. Cane sugar
2. 1N HCl

Procedure

Prepare 20% solution of pure cane sugar by dissolving 20 gms cane sugar in 100 ml distilled water. Pipette out 25 ml of the cane sugar solution in a dry conical flask. Take 25ml 1N HCl solution in a clean dry beaker. Keep both containers i.e. the conical flask (containing sugar solution) and the beaker (containing 1N HCl solution) in a water bath so that the two solutions attain the same temperature.

Mix the cane sugar solution and 1N HCl solution, shake the mixture and immediately fill it in the observation tube; and measure the angle of rotation. Consider this reading as the zero reading.

Take readings for angle or rotation at an interval of 10 minutes.

Take six readings other than the zero reading. Tabulate your readings as shown below

Observation Table

No.	Time in Mins	Rotation angle 'r'	Inversion $r_0 - r_t$	$r_t - r_\infty$	$\log (r_t - r_\infty)$	k
1	0					
2	10					
3	20					
4	30					
5	40					
6	50					
7	60					

Ask the Examiner for infinity reading i.e. r_∞

- ✓ Give the equation to calculate the rate of the reaction
- ✓ Plot graph of $(r_0 - r_t)$ against t
- ✓ Plot graph of $\log(r_t - r_\infty)$ against t
- ✓ Deduce the order of the inversion of cane sugar from the graphs

Result

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

1. The rate constant for the inversion of cane sugar, $k =$ _____
2. From the graphs the order of the reaction is _____

BHAKTA KAVI NARSINH MEHTA UNIVERSITY
BSc Semester V
Chemistry Practical Examination

Centre : _____

Place : _____

The purpose of the practical examination, the division of the batches and distribution of the practical work, candidates are divided into the following groups:

DATES	SEAT NOS. OF BATCH		GROUPS	
	FROM	TO	'A' GROUP	'B' GROUP

DAY	Session	'A' GROUP	'B' GROUP
FIRST DAY	10.00 AM, Session-1	Organic Separation	Physico Chemical
	Session-2	Physico Chemical	Organic Separation
SECOND DAY	Session-3	Volumetric Analysis & Viva	

Note:

The Candidates are informed to be present at the examination centre at least 15 mins before the commencement of the Examination. The following are to be brought by the candidate at the time of the Examinations:

BHAKTA KAVI NARSINH MEHTA UNIVERSITY

BSc Semester V

Chemistry Practical Examination

- Certified Journal
- Hall ticket
- Fee Receipt
- College Identity Card
- Apron / Lab Coat
- Calculator
- Match Box
- Small Knife
- Cloth Duster

Senior Examiner

TYBSc Chemistry Practical Examination

March / April / Oct / Nov 20

